
ERGOSURG

**EM-
Shaverbladetracker
40820123EU**

**EM-shaver blade
tracker 40820123EU**

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1 Intended use

The EM shaver blade tracker allows you to determine the position and orientation of compatible shaver attachments and sinus burrs (see section 'Compatibility'). The EM shaver blade tracker may only be used in conjunction with the electromagnetic navigation system EMU. Failure to observe and follow this instruction manual and the instruction manuals of products used in combination can result in injury to patients, users and third parties as well as damage to the product. Please read all relevant instruction manuals carefully and always follow the instructions given precisely. Check the functioning of the products used in combination. Incorrect application of medical instruments poses a risk of injury for patients. Users of medical instruments must have an appropriate medical qualification and be acquainted with the application.

2 Handling

It is recommended that the suitability of the products for the intended procedure be checked prior to use.

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3 Compatibility

The EM shaver blade tracker is compatible with the following

shaver attachments:

- KARL STORZ shaver attachment (article number 41201 KK)
- KARL STORZ shaver attachment (article number 41204 KKB)
- KARL STORZ sinus burr (article number 41305 D)
- KARL STORZ sinus burr (article number 41305 DW)
- KARL STORZ sinus burr (article number 41303 DT)

NOTE: The EM shaver blade tracker must only be used in combination with the attachments listed in this section.

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4 Safety instructions

Failure to observe and follow this instruction manual and the instruction manuals of products used in combination can result in injury to patients, users and third parties as well as damage to the product. Please read all relevant instruction manuals carefully and always follow the instructions given precisely. Check the functioning of the products used in combination. Incorrect application of medical instruments poses a risk of injury for patients. Users of medical instruments must have an appropriate medical qualification and be acquainted with the application.

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5 Inspection

Incorrectly assembled and damaged instruments can lead to injuries to the patient. Instruments and all of the accessories used in combination with them must be checked immediately before and after every use to ensure that they are complete, free from damage, and in full working order and have no unintentional rough surfaces, sharp corners, burred edges or projecting parts. When connecting the instrument, check that its positioning is stable by pulling on it lightly. Damaged, bent, or incomplete instruments or accessories must not be used. Care must be taken not to leave missing or broken-off components inside the patient.

6 Operating and storage conditions

Always handle the shaver blade tracker with care to avoid damaging it:

- The cable must not be bent.
- Always grasp the plug when removing the shaver blade tracker from the system. Do not pull on the cable.
- Do not subject the cable to excess tension and, in particular, do not stretch the cable.
- Only connect the shaver blade tracker to the designated system sockets.
- The shaver blade tracker must not be exposed to a strong magnetic field, in particular MR devices.
- The shaver blade tracker must not be placed directly on the field generator.
- During operation, the shaver blade tracker's cable must not be placed in the direct vicinity (<30 mm) of the field generator's cable.
- Overloading the instrument by exerting too much force may cause the medical device to break, bend, and malfunction, and consequently injure the patient or user. Do not overload the instruments. Do not bend bent instruments back to their original position.
- Ensure the sterility of the instruments and surroundings.

NOTE: The EM shaver blade tracker must not be used if it is damaged.

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2

7 Connection of the instruments

Remove the EM instruments and the EM patient reference from the sterile packaging. Connect the patient reference to the *applied parts' socket marked 1 on the front of the ETU*. Continue to connect the other instruments to be navigated to the *applied parts' sockets marked '2' to '4' on the front of the ETU*. A status light above the socket shows the connection status:

- Off: no applied part connected/recognized
- Orange: applied part connected but not initialized
- Green: applied part initialized. To initialize the applied parts, switch to the navigation system in the treatment mode.

NOTE: Ensure that at least one patient reference and one probe are connected to the system, otherwise navigation is not possible.

8 Assembly

CAUTION: The tracker can only be mounted in the mounting direction shown. Forcibly mounting against this direction can result in damage to the tracker and/or the shaver blades.

1. Remove the EM shaver blade tracker from the sterile packaging.
2. Slide the shaver blade tracker over the distal end of the shaver attachment and secure it by turning it in the direction of the orientation pin of the set wheel. The safety latch of the tracker must click noticeably into place on the orientation pin.

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9 Disassembly

CAUTION: The tracker can only be removed in the disassembly direction shown. Forcibly removing against this direction can result in damage to the tracker and/or the shaver blades.

1. The EM shaver blade tracker may only be reprocessed when disassembled.
2. Firstly, turn the shaver blade tracker to loosen it from the set wheel's orientation pin.
3. Then pull the shaver blade tracker off the shaver attachment in the direction of the distal end.

10 Reprocessing

The electromagnetic navigation instruments comprise three permanently connected components: the applied part, the cable and the plug.

10.1 Cleaning, disinfection and sterilization before use

These instruments are not sterile when delivered. The use of non-sterile instruments poses a risk of infection for patients, users and third parties. Inspect instruments for visible contamination. Visible contamination is an indication that reprocessing has not been carried out or has been carried out incorrectly.

Reprocess the instruments before initial use and before and after every subsequent use using validated procedures.

Incorrectly reprocessed instruments expose patients, users and third parties to a risk of infection as well as the danger of an instrument malfunction. Observe these instructions and the accompanying documentation.

Effective cleaning and disinfection is a key prerequisite for effective sterilization.

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10.1.1 Special notes

Within the scope of your responsibility for the sterility of the instruments during use, please ensure

- that only sufficiently device- and product-specific validated methods for cleaning/disinfection and sterilization are employed,
- that the devices used (disinfector, sterilizer) are maintained and checked on a regular basis and that the validated parameters are complied with during every cycle.

Please also note the applicable legal regulations in your country as well as the hygiene regulations in the practice or hospital. This applies in particular for the various provisions regarding effective prion inactivation. Reprocessing may only be carried out by suitably trained staff.

NOTE: The use of chemicals which have not been approved may cause damage to the instruments. Only use approved chemicals for reprocessing – see section 'material stability'.

10.2 Cleaning and disinfection**10.2.1 Basics**

Wherever possible, a machine procedure (disinfector) should be used for cleaning and disinfection. A manual procedure should only be used where a machine procedure is unavailable as such procedures are considerably less effective and reproducible.

In both cases preliminary treatment is required.

10.2.2 Preliminary treatment

Heavy soiling must be removed from the instruments directly after use (within a maximum of 2 hours).

1. For this, remove the EM instruments / reference from the surgical unit (ETU) by pulling out the instrument plugs.
2. Use running water for the preliminary treatment.
3. For the manual removal of soiling, only use a suitable, soft, disinfected brush (e.g., intradental brush, cleaning brush KARL STORZ No. 27652) or a soft, clean cloth which you exclusively use for this purpose; do not use metal brushes or steel wool under any circumstances.

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10.3 Machine cleaning/disinfection**10.3.1 Selection of the disinfectant**

When selecting a disinfectant, it must be ensured

- that the disinfectant offers tested effectiveness (e.g., DGHM or FDA approval or CE marking as per DIN EN ISO 15883/ANSI/AAMI ST15883),
- that, if possible, a tested program for thermal disinfection is used (A_0 value > 3000 or, with older devices, at least 5 min at 90 °C) (with chemical disinfection there is a risk of disinfection residues on the instruments),
- that the program used is suitable for the instruments and includes a sufficient number of rinsing cycles,
- that only water which is sterile/almost sterile (max. 10 germs/ml) and low in endotoxins (max. 0.25 endotoxin units/ml) is used for rinsing (e.g., purified water/highly purified water),
- that the air used for drying is filtered and
- that the disinfectant is maintained and checked on a regular basis.

**10.3.2 Selection of the cleaning agent
system**

When selecting the cleaning agent system to be used, it must be ensured

- that this is, as a basic principle, suitable for the cleaning of instruments made from metals, silicones and plastics,
- that, if thermal disinfection is not used, an additional suitable disinfectant which has been proven to be effective (e.g., VAH/DGHM or FDA/EPA approval/ registration or CE marking) is used and that this is compatible with the cleaning agent used and
- that the chemicals used are compatible with the instruments (see section on 'Material stability').
- Cleaning in an ultrasound bath is not permitted.

The concentrations of the cleaning agent and, if relevant, disinfectant, as specified by the manufacturer, must always be adhered to.

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10.3.3Procedure

1. Place the instruments into the disinfectant and when doing so, ensure that the instruments do not touch.
2. Start the program.
3. When the program has finished, remove the instruments from the disinfectant.
4. Check and pack the instruments as promptly as possible after removal (see section on 'Inspection', 'Maintenance' and 'Packaging') in a clean location, after additional drying if necessary.

**10.3.4Fundamental suitability for
effective machine cleaning and
disinfection.**

The fundamental suitability of the instruments for effective machine cleaning and disinfection was evaluated by an independent, accredited testing laboratory using the disinfectant G 7836 CD (thermal disinfection, Miele & Cie. GmbH & Co., Gütersloh, Germany) and the cleaning agent Sekumatic MultiClean (Ecolab, Düsseldorf, Germany).

In so doing, the procedure described above was taken into consideration.

10.4 Manual cleaning

When selecting the cleaning agents and disinfectants to be used, it must be ensured

- that these are, as a basic principle, suitable for the cleaning and disinfection of instruments made from metals, silicones and plastics,
- that a disinfectant which has been proven to be effective (e.g., VAH/DGHM or FDA/EPA approval or CE marking) is used and that this is compatible with the cleaning agent used and
- that the chemicals used are compatible with the instruments (see section on 'Material stability').

**10.4.1Selection of the cleaning agents
and disinfectants**

The use of combined cleaning agents/disinfectants should be avoided as far as possible. Combined cleaning agents/disinfectants may only be used in cases of very low contamination (no visible soiling).

It is imperative that the concentrations and exposure times of the cleaning agents and disinfectants, as specified by the manufacturer are adhered to, in addition to the guidelines for rinsing.

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Only use freshly prepared solutions and only use water which is sterile/almost sterile (max. 10 germs/ml) and low in endotoxins (max. 0.25 endotoxin units/ml) e.g., purified water/highly purified water. In addition, only use filtered air for drying.

10.4.2 Cleaning procedure

1. Place the instruments in the cleaning bath for the specified exposure time so that these are sufficiently covered (if necessary brush carefully with a suitable, soft, disinfected cleaning brush, e.g., intradental brushes, KARL STORZ No. 27652). When doing so, however, ensure that the instruments do not touch. If applicable, rinse all the instruments' lumina at least five times at the start or at the end of the exposure time with a disposable syringe (minimum volume 10 ml), if necessary also use a suitable disposable cannula.
2. Then remove the instruments from the cleaning bath and rinse these thoroughly with water at least three times. If applicable, rinse all the instruments' lumina at least five times with a disposable syringe (minimum volume 10 ml), if necessary also use a suitable disposable cannula.
3. Inspect the instruments (see section on 'Inspection' and 'Maintenance').

10.4.3 Disinfection procedure

1. Immerse the cleaned and inspected instruments in the disinfection bath for the prescribed exposure time in such a way that they are sufficiently covered. When doing so, however, ensure that the instruments do not touch. If applicable, rinse all the instruments' lumina at least five times at the start or at the end of the exposure time with a disposable syringe (minimum volume 10 ml).
2. Then remove the instruments from the disinfection bath and rinse these thoroughly with water at least five times. If applicable, rinse all the instruments' lumina at least five times with a disposable syringe (minimum volume 10 ml).
3. Use filtered compressed air to blow the interior and exterior of the instruments dry.
4. Pack the instruments as promptly as possible after removal (see section on 'Packaging') in a clean location, after additional drying if necessary.

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10.4.4 Fundamental suitability for effective cleaning and disinfection

The fundamental suitability of the EM instruments for effective manual cleaning and disinfection was evaluated by an independent, accredited testing laboratory using the cleaning agent Cidezyme/Enzol and the disinfectant Cidex OPA (Johnson & Johnson GmbH, Norderstedt, Germany). In so doing, the procedure described above was taken into consideration.

10.5 Inspection

After cleaning or cleaning and disinfection, inspect all the EM instruments for corrosion, damaged surfaces, chipping and soiling and single out any damaged instruments (for details on restrictions as to the number of times the instruments can be used, see the section on 'Reusability').

Instruments which are still soiled must be cleaned and disinfected again.

10.6 Maintenance

Instrument oils must not be used.

10.7 Packaging

Please pack the instruments individually in disposable sterilization packaging (single or double packaging) and/or sterilization containers which comply with the following requirements:

- DIN EN ISO/ANSI AAMI ISO 11607
- Suitability for steam sterilization (temperature-resistant up to at least 138°C (280°F), sufficiently vapor permeable)
- Suitable protection of the instruments or sterilization packaging against mechanical damage
- Regular maintenance in accordance with the manufacturer's specifications (sterilization container)

10.8 Sterilization

10.8.1 Sterilization methods

DANGER: Flash sterilization is, in principle, not permitted. Hot air, radiation, formaldehyde, ethylenoxide and plasma sterilization are also not permitted.

Only the sterilization methods listed below may be used for sterilization; other methods are not permitted.

- Fractionated vacuum or gravity displacement method*.
- Steam sterilizer as per DIN EN 13060 or DIN EN 285/ANSI AAMI ST79

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**or 18 min. (prion inactivation – not applicable in the USA)

10.8.2 Fundamental suitability for effective steam sterilization

The fundamental suitability of the instruments for effective steam sterilization was verified by an independent, accredited testing laboratory using the steam sterilizer HST 6x6x6 (Zirbus Technology GmbH, Bad Grund, Germany) and also using the fractionated vacuum procedure as well as the gravity displacement method.

In so doing, the typical conditions in the hospital and doctor's office, as well as the procedure described above, were all taken into consideration.

11 Storage

Following sterilization, the instruments must be stored in their sterile packaging in a dry and dust-free environment.

12 Material stability

When selecting the cleaning agents and disinfectants, it must be ensured that they do not contain the following components:

- Organic, mineral or oxidizing acids (minimum permissible pH value 5.5)
- Strong alkalis (maximum permissible pH value 10, mild alkaline cleaning agents are recommended)
- Organic solvents (e.g. alcohols, ethers, ketones, benzines)
- Oxidizers (e.g., hydrogen peroxides)
- Halogens (chlorine, iodine, bromine)
- Aromatic hydrocarbons/alkyl halides

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- Heavy metal salts
- Oils
- Morpholine
- Never clean the instruments with metal brushes or steel wool.
- The instruments must never be exposed to temperatures in excess of 142°C (288°F)!

13 Directive compliance

This medical device bears the CE mark in accordance with the Medical Device Directive (MDR) 2017/745. A code number after the CE mark indicates the responsible notified body.

14 Reusability

ERGOSURG guarantees a life cycle of 30 reprocessings/uses. If handled with due care and assuming they are neither damaged nor soiled, these instruments can be used again if they still function correctly; any other reuse or the use of damaged and/or soiled instruments is the responsibility of the user. Disregard results in the exclusion of any liability.

15 Disposal

For disposal, no special measures are necessary. National laws and regulations must be observed.

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